

FINAL REPORT OF LABORATORY EXAMINATION
 Humboldtstraße 2, D 70806 Kornwestheim, Germany
 + 49 (0)7154 17 50585

idxbioanalytics-europe@idexx.com www.idexbioanalytics.eu/

IDEXX BioAnalytics Case # 46097-2021

Received Date: 3/12/2021

Completed: 3/17/2021

SEROLOGY DETAILS

	Baseline	1M_TAZ D	2M_TAZ D	3M_TAZ D	4M_TAZ D	5M_TAZ D	6M_TAZ D	7M_TAZ D	8M_TAZ D	9M_TAZ D	10M_TAZ D
Ectromelia											
Ectromelia purified virus	MFI > 2.000	-	-	-	-	-	-	-	-	-	-
EDIM											
EDIM purified virus	MFI > 4.000	-	-	-	-	-	-	-	-	-	-
EDIM VP6 recombinant	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
LCMV											
LCMV recombinant	MFI > 3.000	-	-	-	-	-	-	-	-	-	-
Mycoplasma pulmonis											
<i>M. pulmonis</i> purified bacteria	MFI > 4.250	-	-	-	-	-	-	-	-	-	-
MAV1											
MAV1 purified virus	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
MAV2											
MAV2 purified virus	MFI > 4.500	-	-	-	-	-	-	-	-	-	-
MHV											
MHV purified virus	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
MHV NC recombinant	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
MNV											
MNV purified virus	MFI > 3.000	-	-	-	-	-	-	-	-	-	-
MNV VP1 recombinant	MFI > 4.000	-	-	-	-	-	-	-	-	-	-
MPV											
MPV VP2 recombinant	MFI > 2.000	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
MVM											
MVM VP2 recombinant	MFI > 2.000	-	-	-	-	-	-	-	-	-	-
NS1 ¹	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
PVM											

If you have questions, please call +49 (0)7154 17 50585

	Baseline	1M_TAZ D	2M_TAZ D	3M_TAZ D	4M_TAZ D	5M_TAZ D	6M_TAZ D	7M_TAZ D	8M_TAZ D	9M_TAZ D	10M_TAZ D
PVM purified virus	MFI > 2.500	-	-	-	-	-	-	-	-	-	-
REO3											
REO3 purified virus	MFI > 4.900	-	-	-	-	-	-	-	-	-	-
TMEV											
TMEV purified virus	MFI > 2.250	-	-	-	-	-	-	-	-	-	-
RTV purified virus	MFI > 2.750	-	-	-	-	-	-	-	-	-	-
Sendai											
Sendai purified virus	MFI > 3.750	-	-	-	-	-	-	-	-	-	-
<i>Clostridium piliforme</i>											
<i>Clostridium piliforme</i> purified bacteria	MFI > 4.950	-	-	-	-	-	-	-	-	-	-

NS1¹: NS1 protein is highly conserved among rodent parvoviruses and thus serves as a generic assay for parvovirus seroconversion.

Legend: + = positive - = negative blank = test not performed EQ = equivocal HE = hemolysis precluded testing I = insufficient W = weak positive WB = Western Blot confirmatory analysis pending NS = non-specific reactivity N = normal IgG L = less than normal IgG

Positive MFI results are reported as "+" followed by a number from 1 to 33 in thousands rounded off to the nearest thousand.